



# Anti-Saporin (full length) polyclonal antibody (DPATB-H82815)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-Saporin Polyclonal Antibody
<b>Target</b>	Saporin
<b>Immunogen</b>	Full length protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	Saponaria officinalis
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, IHC-Fr
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Saporin is obtained from the seeds of the Soapwort plant (*Saponaria officinalis*), a plant that grows wildly in Britain and other parts of Europe. Saporin is a plant enzyme with N-glycosidase activity that depurinates a specific nucleotide in the ribosomal

---

**Keywords**

Saporin,

---