



Anti-S. aureus Polyclonal antibody (DPATB-H82817)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Staphylococcus aureus Polyclonal Antibody
Target	S. aureus
Immunogen	Staphylococcus aureus strain ATCC 27660 (mixture of intact organisms and lysed cells).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. aureus
Purification	IgG purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-Fr
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Staphylococcus aureus are gram-positive bacteria that, although part of the normal flora of humans (nasal passages, skin and mucous membranes), have the potential to become opportunistic pathogens. S. aureus is a frequent cause of a wide range of infectio

Keywords

S. aureus,
