



Anti-S. aureus Polyclonal antibody (DPATB-H82818)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Staphylococcus aureus Polyclonal Antibody
Target	S. aureus
Immunogen	Native Staphylococcus aureus, ATCC 27660.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. aureus
Purification	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	ICC/IF
Format	Liquid
Size	1 ml
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Staphylococcus aureus are gram-positive bacteria that, although part of the normal flora of humans (nasal passages, skin and mucous membranes), have the potential to become opportunistic pathogens. S. aureus is a frequent cause of a wide range of infectio

Keywords

S. aureus,
