



Anti-S. aureus Polyclonal antibody (DPATB-H82819)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|--|
| Product Overview | Rabbit Anti-Staphylococcus aureus Polyclonal Antibody |
| Target | S. aureus |
| Immunogen | Soluble and structural antigens of Staphylococcus aureus whole bacterium; ATCC #27660. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | S. aureus |
| Purification | Whole antiserum |
| Conjugate | Unconjugated |
| Applications | ICC/IF, IHC-Fr |
| Format | Liquid |
| Size | 500 µl |
| Buffer | PBS, 10mg/ml BSA |
| Preservative | 0.1% Sodium Azide |
| Storage | Store at +4°C. |

BACKGROUND

Introduction

Staphylococcus aureus are gram-positive bacteria that, although part of the normal flora of humans (nasal passages, skin and mucous membranes), have the potential to become opportunistic pathogens. S. aureus is a frequent cause of a wide range of infectio

Keywords

S. aureus,
