



Anti-S. aureus Polyclonal antibody (DPATB-H82820)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Staphylococcus aureus Polyclonal Antibody
Target	S. aureus
Immunogen	Soluble and structural antigens of Staphylococcus aureus whole bacterium; ATCC #27660.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. aureus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-Fr
Format	Liquid
Size	500 µl
Buffer	PBS, 10mg/ml BSA
Preservative	None
Storage	Store at +4°C.

BACKGROUND

Introduction

Staphylococcus aureus are gram-positive bacteria that, although part of the normal flora of humans (nasal passages, skin and mucous membranes), have the potential to become opportunistic pathogens. S. aureus is a frequent cause of a wide range of infectio

Keywords

S. aureus,
