



# Anti-Staphylococcus type T8 T14/49 T25 TImp.19 protein Group A Polyclonal antibody (DPATB-H82822)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Streptococcus Group A T8+T14/49+T25+TImp.19 protein Polyclonal Antibody
<b>Target</b>	Staphylococcus type T8 T14/49 T25 TImp.19 protein Group A
<b>Immunogen</b>	Tissue / cell preparation of Streptococcus Group A T8+T14/49+T25+TImp.19 protein serotype (heat inactivated)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Staphylococcus
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Agg
<b>Positive Control</b>	Antigenic suspension
<b>Format</b>	Liquid
<b>Size</b>	2 ml
<b>Buffer</b>	Whole serum
<b>Preservative</b>	0.1% Sodium Azide

**Storage**

Store at +4°C. Do not freeze.

---

## BACKGROUND

**Introduction**

Streptococcus is a genus of spherical Gram-positive bacteria belonging to the phylum Firmicutes and the lactic acid bacteria group. Hemolytic Streptococcus are grouped into 23 groups from A to V excluding I and J, according to the serological differences of

---

**Keywords**

GAS, Streptococcus pyogenes,

---