



Anti-Staphylococcus type T9 T18 T22 T23 protein Group A Polyclonal antibody (DPATB-H82823)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus Group A T9+T18+T22+T23 protein Polyclonal Antibody
Target	Staphylococcus type T9 T18 T22 T23 protein Group A
Immunogen	Tissue / cell preparation of Streptococcus Group A T9+T18+T22+T23 protein serotypes (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Staphylococcus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	0.1% Sodium Azide

Storage

Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Streptococcus is a genus of spherical Gram-positive bacteria belonging to the phylum Firmicutes and the lactic acid bacteria group. Hemolytic Streptococcus are grouped into 35 groups from A to V excluding I and J, according to the serological differences of??
Keywords	GAS, Streptococcus pyogenes,
