



Anti-CRPR (aa 1611-1710) polyclonal antibody (DPATB-H82825)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-CRPR Polyclonal Antibody
Target	CRPR
Immunogen	Recombinant fusion protein (tagged): TKYGS LDSGL NRAGIAAPGT NKHTVEG SNP IFNEAIKTPD LDAISEGSND SDLIGIEDLP HFGNVFMDPE VNEKANGYPE VANHNNNF AF NPTPFSPEFV , corresponding to amino acids 1611-1710 of Manduca sexta (tobacco hornworm), CRPR TKYGS LDSGLNRAGIAAPGTNKHTV
Isotype	IgG
Source/Host	Mouse
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Cadherin-related protein receptor BT-R1 (CRPR), belongs to the Cry1A family of Bacillus thuringiensis insecticid
Keywords	Cadherin related protein receptor BT R1;