



Anti-Saporin (full length) polyclonal antibody (DPATB-H82833)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Saporin Polyclonal Antibody
Target	Saporin
Immunogen	Full length protein.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P, IHC-Fr
Format	Liquid
Size	200 µl
Buffer	PBS, Whole serum
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Saporin is obtained from the seeds of the Soapwort plant (*Saponaria officinalis*), a plant that grows wild in Britain and other parts of Europe. Saporin is a plant enzyme with N-glycosidase activity that depurinates a specific nucleotide in the ribosomal

Keywords

Saporin,
