



Anti-DEBS1 (aa 1460-1559) polyclonal antibody (DPATB-H82836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-DEBS1 Polyclonal Antibody
Target	DEBS1
Immunogen	Recombinant fragment: RPTGAGEPAR LDGTWLVAKY AGTADETSTA AREALESAGA RVRELVDAR CGRDELAERL RSVGEVAGVL SLLAVDEAEP EEAPLALASL ADTSLVQAM , corresponding to amino acids 1460-1559 of DEBS1 (Saccharopolyspora erythraea) RPTGAGEPARLDGTWLVAKYAGTADETSTAAREALESAGA RV
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Saccharopolyspora erythraea
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

This antibody was raised by a genetic immunization technique. Genetic immunization can be used to generate antibodies by directly delivering antigen-coding DNA into the animal, rather than injecting a protein or peptide (Tang et al. PubMed: 1545867; Chamb

Keywords

EryAI; EryAI Erythromycin polyketide synthase modules 1 and 2;
