



Anti-S. boydii type 15 Polyclonal antibody (DPATB-H82837)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Shigella boydii serotype 15 Polyclonal Antibody
Target	S. boydii type 15
Immunogen	Tissue / cell preparation of Shigella boydii serotype 15 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. boydii
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Shigella is a nonmotile gram-negative bacillus belonging to the Enterobacteriaceae family. The organisms taken in orally invade the epithelial cell of the large intestine to cause shigellosis with (mucous and bloody) diarrhea, fever and abdominal pain as

Keywords

S.boydii,
