



Anti-S. boydii type 4 Polyclonal antibody (DPATB-H82838)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Shigella boydii serotype5 Polyclonal Antibody
Target	S. boydii type 4
Immunogen	Tissue / cell preparation of Shigella boydii serotype5(heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. boydii
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Shigella is a non motile gram negative bacillus belonging to the Enterobacteriaceae family. The organisms invade the epithelia of the large intestine to cause shigellosis. Shigella is classified into 5 species based on the serological properties of the som

Keywords

S.boydii serotype 4,
