



Anti-Streptokinase (native) polyclonal antibody (DPATB-H82843)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptokinase Polyclonal Antibody
Target	Streptokinase
Immunogen	Full length native protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	500 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Streptokinase is a metallo-enzyme produced by various strains of Streptococcus that can liquify
--------------	---

blood clots by converting plasminogen to plasmin. Streptokinase is not a protease but it activates plasminogen by complexing with it. It is thought to prevent

Keywords

ska,skc,skg,
