



Anti-Staphylococcus Group A Polyclonal antibody (DPATB-H82844)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus group A Polyclonal Antibody
Target	Staphylococcus Group A
Immunogen	Type specific carbohydrate of group A Streptococcus
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Staphylococcus
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ICC/IF, ELISA
Format	Liquid
Size	500 µl
Buffer	PBS, 10mg/ml BSA
Preservative	0.1% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction

The genus Streptococcus is comprised of a wide variety of both pathogenic and commensal gram positive bacteria which are found to inhabit a wide range of hosts, including humans, horses, pigs and cows. Within the host, streptococci are often found to colo

Keywords

Group A strep, Group A streptococci, Group A streptococcus,
