



Anti-S. pneumoniae type 1, 2, 3, 4, 5 Polyclonal antibody (DPATB-H82845)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus pneumoniae type 1+2+3+4+5 Polyclonal Antibody
Target	S. pneumoniae type 1, 2, 3, 4, 5
Immunogen	Tissue / cell preparation of Streptococcus pneumoniae type 1+2+3+4+5 (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pneumoniae
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction Bronchial pneumonia is most prevalent in infants, young children and aged adults. It is caused by various bacteria, including *Streptococcus pneumoniae*. *Streptococcus pneumoniae* are Gram positive, lancet shaped cocci (elongated cocci with a slightly pointed

Keywords *S. pneumoniae*, *Streptococcus pneumoniae*,
