



Anti-S. pneumoniae type 15 Polyclonal antibody (DPATB-H82846)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus pneumoniae type 15 Polyclonal Antibody
Target	S. pneumoniae type 15
Immunogen	Tissue / cell preparation of Streptococcus pneumoniae type 15 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pneumoniae
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Bronchial pneumonia is most prevalent in infants, young children and aged adults. It is caused by various bacteria, including <i>Streptococcus pneumoniae</i> . <i>Streptococcus pneumoniae</i> are Gram positive, lancet shaped cocci (elongated cocci with a slightly pointed
Keywords	<i>S. pneumoniae</i> type 15,