



Anti-S. pneumoniae type 20 Polyclonal antibody (DPATB-H82848)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus pneumoniae type 20 Polyclonal Antibody
Target	S. pneumoniae type 20
Immunogen	Tissue / cell preparation of Streptococcus pneumoniae type 20 (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pneumoniae
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Streptococcus pneumoniae, or pneumococcus, is a Gram positive, alpha hemolytic diplococcus aerotolerant anaerobe and a member of the genus *Streptococcus*. A significant human pathogenic bacterium, *S. pneumoniae* was recognized as a major cause of pneumonia i

Keywords

S. pneumoniae type 20,