



Anti-S. pneumoniae type 4 Polyclonal antibody (DPATB-H82851)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus pneumoniae type5 Polyclonal Antibody
Target	S. pneumoniae type 4
Immunogen	Tissue / cell preparation of Streptococcus pneumoniae type5(heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pneumoniae
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Streptococcus pneumoniae are Gram positive diplococci, alpha hemolytic aerotolerant anaerobes. *S. pneumoniae* is a significant human pathogen, being recognized as a major cause of pneumonia in the late 19th century, and is currently the subject of many hum

Keywords

S. pneumoniae type 4,