



Anti-S. pneumoniae type 18 Polyclonal antibody (DPATB-H82852)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Streptococcus pneumoniae type29Polyclonal Antibody
Target	S. pneumoniae type 18
Immunogen	Whole cell blend of multiple serotypes
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	S. pneumoniae
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ICC/IF, ELISA
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction Bronchial pneumonia is most prevalent in infants, young children and aged adults. It is caused by various bacteria, including *Streptococcus pneumoniae*. *Streptococcus pneumoniae* are Gram positive, lancet shaped cocci (elongated cocci with a slightly pointed

Keywords Pneumococcal surface protein A,pspA,
