



# Anti- Streptokinase Polyclonal antibody (DPATB-H82853)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-Streptokinase Polyclonal Antibody
<b>Target</b>	Streptokinase
<b>Immunogen</b>	Native, purified Streptokinase.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	Streptococcus.
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB, FuncS
<b>Format</b>	Liquid
<b>Size</b>	1 ml
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Streptokinase is a metallo-enzyme produced by various strains of Streptococcus that can liquify blood clots by converting plasminogen to plasmin. Streptokinase is not a protease but it activates plasminogen by complexing with it. It is thought to prevent

---

**Keywords**

ska,skc,skg,

---