



# Anti-T. pallidum Polyclonal antibody (DPATB-H82861)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Treponema pallidum Polyclonal Antibody
<b>Target</b>	T. pallidum
<b>Immunogen</b>	Highly purified Treponema pallidum.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	T. pallidum
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, ICC/IF
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	10mg/ml BSA, 0.01M PBS
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Treponema pallidum is the causative agent of syphilis. It is a gram negative spirochete, a helical to sinusoidal bacterium with outer and cytoplasmic membranes, a thin peptidoglycan layer, and periplasmic flagella.

---

**Keywords**

T. pallidum,

---