



Anti-T. pallidum Polyclonal antibody (DPATB-H82862)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Treponema pallidum Polyclonal Antibody
Target	T. pallidum
Immunogen	Tissue / cell preparation-highly purified Treponema pallidum
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	T. pallidum
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Size	500 µl
Buffer	PBS, 10mg/ml BSA
Preservative	None
Storage	Store at +4°C.

BACKGROUND

Introduction

Treponema pallidum is the causative agent of syphilis. It is a gram negative spirochete, a helical to sinusoidal bacterium with outer and cytoplasmic membranes, a thin peptidoglycan layer, and periplasmic flagella.

Keywords

T. pallidum,
