



Anti-V. parahaemolyticus K19 Polyclonal antibody (DPATB-H82864)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Vibrio parahaemolyticus K19 Polyclonal Antibody
Target	V. parahaemolyticus K19
Immunogen	Tissue / cell preparation of Vibrio parahaemolyticus K19 (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	V. parahaemolyticus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Vibrio parahaemolyticus is a non spore forming gram negative bacillus with monopolar flagella. This bacterium is frequently isolated from estuarine and marine environments and is a leading cause of seafood associated gastroenteritis. Serological types are

Keywords

V. parahaemolyticus K19,
