



Anti-*V. parahaemolyticus* K48 Polyclonal antibody (DPATB-H82871)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Vibrio parahaemolyticus K48 Polyclonal Antibody
Target	<i>V. parahaemolyticus</i> K48
Immunogen	Tissue / cell preparation of <i>Vibrio parahaemolyticus</i> K48
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>V. parahaemolyticus</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Vibrio parahaemolyticus is a curved, rod shaped, Gram negative bacterium found in brackish saltwater, which, when ingested, causes gastrointestinal illness in humans. V. parahaemolyticus bacterium is in the same family as those that cause cholera. V. para

Keywords

V. parahaemolyticus K48,