



# Anti-V. parahaemolyticus K49 Polyclonal antibody (DPATB-H82872)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Rabbit Anti-Vibrio parahaemolyticus K49 Polyclonal Antibody
Target	V. parahaemolyticus K49
Immunogen	Tissue / cell preparation of Vibrio parahaemolyticus K49 (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	V. parahaemolyticus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	0.1% Sodium Azide
Storage	Store at +4°C. Do not freeze.

## BACKGROUND

**Introduction**

Vibrio parahaemolyticus is a curved, rod shaped, Gram negative bacterium found in brackish saltwater, which, when ingested, causes gastrointestinal illness in humans. V. parahaemolyticus bacterium is in the same family as those that cause cholera. V. para

**Keywords**

V. parahaemolyticus K49,