



Sheep anti-JWH-018 Polyclonal antibody (DPATB-H82892)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-JWH-018 Polyclonal Antibody
Target	JWH-018
Immunogen	JWH-018 conjugated to Bovine thyroglobulin.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Bovine
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	pH: 7.20Preservative: 0.09% Sodium azide98% Phosphate Buffer, 0.88% Sodium chloride
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

JWH-018 (1-pentyl-3-(1-naphthoyl)indole) is an analgesic chemical from the naphthoylindole family, which acts as a full agonist at both the CB1 and CB2 cannabinoid receptors, with some selectivity for CB2. It produces effects in animals similar to those o

Keywords

1-pentyl-3-(1-naphthoyl)indole;
