



Anti-H. pylori Polyclonal antibody (DPATB-H82902)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Helicobacter pylori, prediluted Polyclonal Antibody
Target	H. pylori
Immunogen	This antibody was developed using Helicobacter pylori strain CH-20 429.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	H. pylori
Conjugate	Unconjugated
Applications	IHC-P, IHC-Fr
Positive Control	Helicobacter pylori infected Human cells or tissue.
Format	Prediluted
Buffer	99% PBS, 1% BSA
Preservative	0.1% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	The spiral shaped bacterium?Helicobacter pylori?is strongly associated with inflammation of the stomach and is also implicated in the development of gastric malignancy.?H. pylori?is known to
---------------------	--

cause peptic ulcers and chronic gastritis in human. It is associ

Keywords

H. pylori;
