



Magic™ Anti-Gentamycin polyclonal antibody (DPATB-H82908)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Gentamycin Polyclonal Antibody
Target	Gentamycin
Immunogen	Gentamycin Albumin conjugate.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	IgY fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	500 µg
Buffer	Tris. pH 7.4
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.

Keywords

Gentamicin;
