



Magic™ Chicken Anti-Protein G polyclonal antibody (DPATB-H82918U)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein G
Isotype	IgY
Source/Host	Chicken
Species Reactivity	N/A
Purification	Affinity purified.
Conjugate	Unconjugated
Applications	ELISA, WB Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Lot specific
Size	1 mg
Buffer	PBS with 0.09% sodium azide
Preservative	0.09% sodium azide
Storage	Liquid solution stored at 2-8°C is stable for a minimum of 1 year.
Ship	Wet ice