



Anti-HIV type 1 Polyclonal antibody (DPATB-H82945)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-HIV1 Polyclonal Antibody
Target	HIV type 1
Immunogen	Human isolate, confirmed.
Isotype	IgG
Source/Host	Goat
Species Reactivity	HIV
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ICC, IF, WB
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2.No stabilizing proteins have been added.
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

HIV is a highly variable virus which mutates very readily. This means there are many different strains of HIV, even within the body of a single infected person. The strains of HIV1 can be classified into three groups : the "major" group M, the "outlier" g

Keywords

HIV 1; HTLV3; Human Immunodeficiency Virus 1;
