



# Anti-FeLV P27 Polyclonal antibody (DPATB-H82956)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti-FeLV p27 Polyclonal Antibody
<b>Target</b>	FeLV P27
<b>Immunogen</b>	Purified FeLV p27.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	FeLV
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Integral membrane protein.
<b>Positive Control</b>	Purified FeLV p27 or whole virus to FeLV
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

The feline leukemia virus (FeLV) is the causative agent of the most important fatal infectious disease complex of domestic cats today. It is an RNA virus belonging to the family Retroviridae.

---

**Keywords**

feline leukemia provirus; Feline leukemia virus p27; Feline leukemia virus subgroup c receptor related protein 1; Feline sarcoma virus; flvcr protein; flvcr1; fps; hflvcr; Oncogene fes;

---