



# Anti-HSV type 2 Polyclonal antibody (DPATB-H82957)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-HSV2 Polyclonal Antibody
<b>Target</b>	HSV type 2
<b>Immunogen</b>	HSV Type3Strain G (Human)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	HSV
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, Conjugation, IM, WB, ICC, IF
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	0.01M PBS, pH 7.2
<b>Preservative</b>	0.1% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Herpes simplex type3(HSV2) belongs to a family that includes HSV1, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus. HSV1 and HSV4 are extremely difficult to distinguish from each other. These viruses have a DNA genome, an icosahedral p??

---

**Keywords**

Herpes Simplex Virus Type 2;

---