



Anti-HIV type 1 Polyclonal antibody (DPATB-H82963)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-HIV1 Polyclonal Antibody
Target	HIV type 1
Immunogen	Human isolate, confirmed
Isotype	IgG
Source/Host	Goat
Species Reactivity	HIV
Purification	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	ICC, IF
Format	Liquid
Size	250 µl
Buffer	10mM PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

HIV is a highly variable virus which mutates very readily. This means there are many different strains of HIV, even within the body of a single infected person. The strains of HIV1 can be classified into three groups : the "major" group M, the "outlier" g

Keywords

HIV 1; HTLV3; Human Immunodeficiency Virus 1;
