



Anti-Prostaglandin E2 polyclonal antibody (DPATB-H82964)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Prostaglandin E2 Polyclonal Antibody
Target	Prostaglandin E2
Immunogen	Highly pure Prostaglandin E2 conjugated to bovine thyroglobulin
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, IHC-P
Positive Control	recombinant PGE2
Format	Liquid
Size	100 µl
Buffer	Contains Sodium Azide as preservative.
Preservative	See individual product datasheet
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Prostaglandin F2α (PGF2α in prostanoid nomenclature), pharmaceutically termed dinoprost (INN), is a naturally occurring prostaglandin used in medicine to induce labor and as an abortifacient. In domestic mammals, it is produced by the uterus when stimulate
Keywords	PG E2; PGE 2; PGE2; ProstaglandinE2;