



# Anti-*P. pastoris* Polyclonal antibody (DPATB-H82972)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Goat Anti- <i>P.pastoris</i> Polyclonal Antibody
<b>Target</b>	<i>P. pastoris</i>
<b>Immunogen</b>	Tissue / cell preparation: Solubilized <i>Pichia pastoris</i> yeast cells.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Goat
<b>Species Reactivity</b>	<i>P. pastoris</i>
<b>Purification</b>	Protein G purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Format</b>	Liquid
<b>Size</b>	100 µg
<b>Buffer</b>	PBS, Sodium citrate
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

*Pichia pastoris*, a methylotrophic yeast, has been used as a heterologous protein expression system. It combines many of the benefits of *E. coli* expression with the advantages of expression in an eukaryotic system such as protein processing, folding, and p

---

**Keywords**

*Pichia pastoris*;

---