



# Anti-Gentamycin polyclonal antibody (DPATB-H82986)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	Goat Anti-Gentamycin Polyclonal Antibody
Target	Gentamycin
Immunogen	Chemical / Small Molecule: Gentamicin.
Isotype	IgG
Source/Host	Goat
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	RIA
Format	Liquid
Size	250 µl
Buffer	Whole serum.
Preservative	0.1% Sodium Azide
Storage	Add glycerol to a final volume of 50%, aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.

---

**Keywords**

Gentamicin;

---