



# Anti-Methaqualone polyclonal antibody (DPATB-H82994)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-Methaqualone Polyclonal Antibody
<b>Target</b>	Methaqualone
<b>Immunogen</b>	Methaqualoneconjugated to BSA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	cELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	20mM PBS, 150mM Sodium chloride, pH 7.2
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Methaqualone is a nonbarbiturate sedative/hypnotic that produces sleep for about 7 to 9 hours. It also produces muscular relaxation, feelings of contentment, and total passivity. It is an addictive drug.

---

**Keywords**

methaqualone;

---