



Anti-Ceftiofur polyclonal antibody (DPATB-H82995)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Ceftiofur Polyclonal Antibody
Target	Ceftiofur
Immunogen	Ceftiofur conjugated to BSA.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Ceftiofur, a semisynthetic cephalosporin, is a broad-spectrum antibiotic against both Gram-positive and Gram-negative bacteria including beta-lactamase-producing bacterial strains and anaerobes. Its antibacterial activity results from the inhibition of mu
Keywords	[6r [6a7b(z)]]8[[[(2 amino5thiazolyl)(methoxyimino)acetyl]amino] 3 [[2 furanylcarbonyl)thio] methyl]9oxo6thia 1 azabicyclo [4.2.0] oct3ene3carboxylic acid;