



Anti-Methadol polyclonal antibody (DPATB-H82998)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Methadol Polyclonal Antibody
Target	Methadol
Immunogen	Methadol-BSA.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at -20°C.

BACKGROUND

Introduction	Methadol is a narcotic analgesic with a long onset and duration of action. It is used mainly in the treatment of narcotic dependence.
Keywords	Methadone;