



# Anti-Methyldihydrotestosterone polyclonal antibody (DPATB-H83009)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Sheep Anti-Methyldihydrotestosterone Polyclonal Antibody
<b>Target</b>	Methyldihydrotestosterone
<b>Immunogen</b>	17alpha-Methyldihydrotestosterone-BSA
<b>Isotype</b>	IgG
<b>Source/Host</b>	Sheep
<b>Species Reactivity</b>	N/A
<b>Purification</b>	IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	cELISA
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	20mM PBS, 150mM Sodium chloride, pH 7.2
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## BACKGROUND

**Introduction**

Methyldihydrotestosterone is an androgenic anabolic steroid.

---

**Keywords**

MDHT;

---