



Anti-KSHV ORF8 Polyclonal antibody (DPATB-H83038)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-KSHV ORF8 Polyclonal Antibody
Target	KSHV ORF8
Immunogen	Synthetic peptide, corresponding to amino acids sequence of KSHV ORF8.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	KSHV
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Positive Control	TPA-induced BCBL-1 cell lysate.
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Kaposi's sarcoma-associated herpesvirus (KSHV) belongs to the gamma-herpesvirus subfamily and has been closely linked to the Kaposi's sarcoma, primary effusion lymphoma (PEL) and multicentric Castleman's disease. The genome of KSHV is 165-170 kb and cont

Keywords

Glycoprotein B-like; HHV 8; Human herpes virus 8; Human herpesvirus 8; Kaposi sarcoma herpes virus; Kaposi sarcoma herpes virus; KSHV; ORF8;
