



Anti-L Kynurenine polyclonal antibody (DPATB-H83039)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit L Kynurenine (conjugated) Polyclonal Antibody
Target	L Kynurenine
Immunogen	Synthetic peptide L Kynurenine conjugated to a protein carrier.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Cellular Localization	Cytoplasmic
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	L-Kynurenine is a metabolite of the amino acid L-tryptophan used in the production of niacin. Kynurenine is synthesized by the enzyme tryptophan dioxygenase, which is made primarily but not exclusively in the liver, and indoleamine 2,3-dioxygenase which i
Keywords	Kynureninase; L kynurenine hydrolase;