



Anti-Quinaldic Acid polyclonal antibody (DPATB-H83040)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Quinaldic Acid (conjugated) Polyclonal Antibody
Target	Quinaldic Acid
Immunogen	Synthetic peptide Quinaldic Acid conjugated to a protein carrier.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction The N-heterocycle quinaldic acid (quinoline 2-carboxylic acid) is utilized by *Alcaligenes* sp. F-2

as sole source of carbon and energy.

Keywords

Quinoline3carboxylic acid;
