



Anti-Pyridoxal polyclonal antibody (DPATB-H83041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Pyridoxal (conjugated) Polyclonal Antibody
Target	Pyridoxal
Immunogen	Synthetic Pyridoxal (reduced form) conjugated to BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IHC-Fr, ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Pyridoxal is one of the three natural organic compounds, along with pyridoxamine and
---------------------	--

pyridoxine, that comprise vitamin B6. Each of them is converted by the liver into a single biologically active form, pyridoxal 5-phosphate.

Keywords

Pyridoxine; Vitamin B6;
