



Anti-D-Pantothenic Acid polyclonal antibody (DPATB-H83042)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit d-Pantothenic Acid (conjugated) Polyclonal Antibody
Target	d-Pantothenic Acid
Immunogen	Synthetic d-Pantothenic Acid conjugated to BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IHC-Fr, ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Pantothenic acid is a member of the B-vitamin family, is an essential nutrient sometimes called

vitamin B5. The biologically active form of the vitamin is termed d-pantothenic acid. It is involved in a number of biological reactions, including energy prod

Keywords

Pantothenic acid; Vitamin B5;
