



Anti-Caproic acid polyclonal antibody (DPATB-H83046)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Caproic acid (conjugated) Polyclonal Antibody
Target	Rat Caproic acid
Immunogen	Synthetic Caproic acid conjugated to BSA
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Preadsorbed antiserum
Conjugate	Unconjugated
Applications	ELISA, IHC-Fr
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Hexanoic acid, also known as caproic acid, is the carboxylic acid derived from hexane with the formula $C_5H_{11}COOH$. It is a colorless oily liquid smelling of cheese. It is a fatty acid found naturally in various animal fats and oils.

Keywords

caproic acid; capronic acid; heaxanoic acid;
