



# Anti-Caprylic acid polyclonal antibody (DPATB-H83047)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rat Caprylic acid (conjugated) Polyclonal Antibody
<b>Target</b>	Rat Caprylic acid
<b>Immunogen</b>	Chemical / Small Molecule, Caprylic acid conjugated to BSA.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rat
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Preadsorbed IgG fraction
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-Fr, ELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	50% Glycerol, ddH2O
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Caprylic acid is the common name for the eight-carbon straight chain fatty acid known by the systematic name octanoic acid. It is found naturally in coconuts and breast milk. It is an oily liquid with a slightly unpleasant rancid taste that is minimally s

---

**Keywords**

Octanoic acid;

---