



# Anti-Histidine polyclonal antibody (DPATB-H83056)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Histidine (conjugated) Polyclonal Antibody
<b>Target</b>	Histidine
<b>Immunogen</b>	Synthetic histidine conjugated to protein carrier.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	N/A
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	50% Glycerol, ddH2O
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Histidine is a non essential amino acid found abundantly in hemoglobin. It has been utilised in the treatment of rheumatoid arthritis, allergic diseases, ulcers and anemia.

---

**Keywords**

H; His;

---