



Anti-L Serine polyclonal antibody (DPATB-H83057)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit L Serine (conjugated) Polyclonal Antibody
Target	L Serine
Immunogen	Synthetic peptide conjugated to BSA by a Glutaraldehyde linker.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Preadsorbed antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

The amino acid L-serine, is one of the so-called non-essential amino acids; it plays a central role in cellular proliferation and is necessary for specific functions in the central nervous system. Defects in the synthesis of L-serine alter the development
