



Anti-Pantoea agglomerans polyclonal antibody (DPATB-H83066)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Pantoea agglomerans (conjugated) Polyclonal Antibody
Target	Pantoea agglomerans
Immunogen	Pantoea agglomerans total proteins.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA, WB
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

P. agglomerans, formerly named *Enterobacter agglomerans*, is a ubiquitous Enterobacteriaceae that is found in plants and in the feces of humans and animals.

Keywords

Enterobacter agglomerans; *Erwinia herbicola*; *Pantoea agglomerans*;
